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July 20th 2021

## This letter:

- Multiple batches of swabs involved in false-positive results for *B. parapertussis*-PCR: 8 batches e-swabs (orange cap) and 2 batches dry swabs (blue cap).

## Notification:

In the Labinf@ct *B. parapertussis* of June 3rd 2021, we reported that the use of e-swabs can lead to false-positive *B. parapertussis*-PCRs. We requested laboratories to share their experiences and findings with us.

The information which we received, and which was shared by the Medical Molecular Microbiologists in the WMDI (working group Molecular Diagnostics in Infectious diseases), led to the identification of *IS1001*-positive swab tips in Copan swabs. These Copan swabs were found both with orange caps (sampling set) and with blue caps (dry swabs).

- E-swabs with orange cap: Six laboratories reported 8 different batch numbers with false-positive PCRs: lot# 21087100; 203038000; 202323100; 202433300; 203193300; 201896300; 202011501; 2031933.
- (dry) swabs with blue caps: lot# 2004732; 202472700.

All false-positive PCRs concern the detection of the *IS1001*-element; Ct values typically are Ct34 or higher. The *IS1002*-element has not been detected up to now.

At least three laboratories have reported positive results to physicians. Notification of *B. parapertussis*-PCRs are only obligatory in cases with pertussis-like symptoms. Additionally,

notifications may be obligatory if multiple cases are found within an institution (article 26, law on public health). For the full list of criteria, see the LCI guideline Pertussis.

Several of the notified cases have been retracted due to false-positive PCRs. In other cases, the possibility of a false-positive PCR cannot be investigated. Consequently, these notifications will not be retracted.

The distributor (MLS), producer (Copan) and the national health authority have been informed. The national health authority will investigate whether patients have incorrectly received treatment. The extent of the problem cannot be determined based on the notifications, as notifications are not mandatory for all positive PCRs.

All e-swab users have received a letter from the distributor MLS and Copan stipulating that Copan has confirmed the possibility of false-positive PCRs on swabs with an orange cap. Moreover, the false-positivity depends on the sensitivity of the used kit and platform. The letter further mentions:

*Copan cannot exclude the presence of traces of nucleic acid materials into the eSwab that may be detected as Bordetella parapertussis under certain testing conditions.*

*As a good clinical practice, diagnostic companies and end users are required to verify the compatibility of the transport system with the specific molecular assay before use.*

[...]


*Copan became aware of this combination during the ongoing investigation.*

As part of the production process, e-swabs are sterilised with ionizing radiation. See: [ESwab-  
Package-Insert.pdf \(copanusa.com\)](#). In line with this, all bacterial cultures are negative, however DNA remains detectable depending on the kit and platform used.

---

**Authors:** Ellen Mascini, Guido Bastiaens, Jacky Flipse, Angelino Tromp vanuit Rijnstate; Erik van Hannen (WMDI); Erik Schaftenaar (AZN); Bettie Voordouw, Titia Kortbeek, George Sips, Michelle Stadlander, Dimphye van Meijeren, Hans van den Kerkhof (RIVM-Cib)

**Sent to:** members of Labinf@ct

**Appendix Figure.** Comparison of notified cases of *Bordetella pertussis* and *B. paraperussis* registered by the National Institute for Public Health and the Environment (RIVM), the Netherlands, 2003–June 2021. A) *Bordetella pertussis*; B) *B. paraperussis*. Annual reported cases are stratified by age groups.