



Algorithm to identify *Candida auris* based on phenotypic laboratory method and initial species identification

PURPOSE

Candida auris is a multidrug-resistant yeast that has been found in multiple countries, including the United States. *C. auris* can cause invasive infections, be passed from person to person, and persist in the environment. Its severity, communicability, and drug resistance makes correctly identifying *C. auris* crucial to treating patients and preventing infections. However, this is challenging because traditional phenotypic methods frequently misidentify *C. auris*. This algorithm details the steps needed to determine the correct *Candida* spp. based on the tests and equipment available in your lab.

TABLE OF CONTENTS – ALGORITHMS BY METHOD

1. Bruker Biotyper MALDI-TOF
2. bioMérieux VITEK MS MALDI-TOF
3. VITEK 2 YST
4. API 20C
5. BD Phoenix
6. MicroScan
7. RapID Yeast Plus
8. Summary of this algorithm in table form

Please note that these algorithms are based on our current knowledge about misidentification of C.auris and may change as we learn new information.

Identification Method


Database/ Software

Initial finding

Determination


Next steps

Bruker Biotyper MALDI-TOF



C. auris confirmed:
Place patient in transmission-based precautions, report to CDC (candidaauris@cdc.gov), and notify state and local health departments.

bioMérieux VITEK MS
MALDI-TOF




***C. auris* confirmed:**

Place patient in transmission-based precautions, report to CDC (candidaauris@cdc.gov), and notify state and local health departments.

***C. auris* possible:**

Further work-up needed to determine if actually *C. auris*. Send isolates to a reference lab, a state public health lab, a regional lab, or CDC for further identification. Place patient in transmission-based precautions and notify state and local health departments and CDC (candidaauris@cdc.gov).


VITEK 2 YST

***C. auris* confirmed:**

Place patient in transmission-based precautions, report to CDC (candidaauris@cdc.gov), and notify state and local health departments.

***C. auris* possible:**

Further work-up needed to determine if actually *C. auris*. Send isolates to a reference lab, a state public health lab, a regional lab, or CDC for further identification. Place patient in transmission-based precautions and notify state and local health departments and CDC (candidaauris@cdc.gov).



C. auris suspected:
Place patient in transmission-based precautions and notify state and local health departments and CDC (candidaauris@cdc.gov). Send any isolates suspected to be *C. auris* to a reference lab, a state public health lab, a regional lab, or CDC for further identification.

Likely not C. auris:
No further *C. auris*-related action required.

Identification Method

Initial finding

Determination

Next steps

BD Phoenix

C. catenulata

C. haemulonii


Candida spp.
not identified

C. auris possible:
Needs further work-up

C. auris possible:
Needs further work-up


C. auris possible:
Needs further work-up

***C. auris* possible:**
Further work-up needed to determine if actually *C. auris*. Send isolates to a reference lab, a state public health lab, a regional lab, or CDC for further identification. Place patient in transmission-based precautions and notify state and local health departments and CDC (candidaauris@cdc.gov).

***C. auris* possible:**

Further work-up needed to determine if actually *C. auris*. Send isolates to a reference lab, a state public health lab, a regional lab, or CDC for further identification. Place patient in transmission-based precautions and notify state and local health departments and CDC (candidauris@cdc.gov).

RapID Yeast Plus



***C. auris* possible:** Further work-up needed to determine if actually *C. auris*. Send isolates to a reference lab, a state public health lab, a regional lab, or CDC for further identification. Place patient in transmission-based precautions and notify state and local health departments and CDC (candidaauris@cdc.gov).

Identification Method	Database/Software, if applicable	<i>C. auris</i> is confirmed if initial identification is <i>C. auris</i> .	<i>C. auris</i> is possible if the following initial identifications are given. Further work-up is needed to determine if the isolate is <i>C. auris</i> .
Bruker Biotyper MALDI-TOF	RUO libraries (Versions 2014 [5627] and more recent)	<i>C. auris</i>	n/a
	CA System library (Version Claim 4)	<i>C. auris</i>	n/a
bioMérieux VITEK MS MALDI-TOF	RUO library (with Saramis Version 4.14 database and Saccharomycetaceae update)	<i>C. auris</i>	<i>C. haemulonii</i> No identification
	IVD library	n/a	<i>C. haemulonii</i> No identification
VITEK 2 YST	Software version 8.01	<i>C. auris</i>	<i>C. haemulonii</i> <i>C. duobushaemulonii</i> <i>Candida</i> spp. not identified
	Older versions	n/a	<i>C. haemulonii</i> <i>C. duobushaemulonii</i> <i>Candida</i> spp. not identified
API 20C		n/a	<i>Rhodotorula glutinis</i> (with characteristic red color present) <i>C. sake</i> <i>Candida</i> spp. not identified
BD Phoenix		n/a	<i>C. catenulata</i> <i>C. haemulonii</i> <i>Candida</i> spp. not identified
MicroScan		n/a	<i>C. lusitaniae</i> * <i>C. guilliermondii</i> * <i>C. parapsilosis</i> * <i>C. famata</i> <i>Candida</i> spp. not identified
RapID Yeast Plus		n/a	<i>C. parapsilosis</i> * <i>Candida</i> spp. not identified

* *C. guilliermondii*, *C. lusitaniae*, and *C. parapsilosis* generally make hyphae or pseudohyphae on cornmeal agar. If hyphae or pseudohyphae are not present on cornmeal agar, the isolate should raise suspicions of being *C. auris* as *C. auris* typically does not make hyphae or pseudohyphae. However, some *C. auris* isolates have formed hyphae or pseudohyphae. Therefore, it would be prudent to consider any *C. guilliermondii*, *C. lusitaniae*, and *C. parapsilosis* isolates identified on MicroScan and any *C. parapsilosis* isolates identified on RapID Yeast Plus as possible *C. auris* isolates and further work-up should be considered.

If *C. auris* is confirmed: Place patient in transmission-based precautions, report to CDC (candidaauris@cdc.gov), and notify state and local health departments.

If *C. auris* is possible: Further work-up is needed to determine if actually *C. auris*. Send isolates to a reference lab, a state public health lab, a regional lab, or CDC for further identification. Place patient in transmission-based precautions and notify state and local health departments and CDC (candidaauris@cdc.gov).