

SUPPLEMENTAL INFORMATION

Trisaminohexyl Isocyanurate, a Urinary Biomarker of HDI Isocyanurate Exposure

Zachary Robbins, Wanda Bodnar, Zhenfa Zhang, Avram Gold, and Leena A. Nylander-French*

Department of Environmental Sciences and Engineering, Gillings School of Global Public Health, The University of North Carolina at Chapel Hill, NC 27599, USA.

Corresponding Author

* Leena A. Nylander-French, Ph.D., CIH, Department of Environmental Sciences and Engineering, Gillings School of Global Public Health, The University of North Carolina at Chapel Hill, NC 27599, USA, Email: leena_french@unc.edu

Table of Contents

Table..... S2
Figures..... S2 – S9

Table S1. Gradient program for the NanoAcquity ultra-performance liquid chromatography system (Waters Corp., Milford, MA) using a Symmetry C18 trapping column, 5 μm , 180 μm \times 20 mm (Waters Corp.) and an Atlantis dC18 analytical column, 3 μm , 100 μm \times 100 mm (Waters Corp.). Mobile phase A consisted of 0.1% formic acid in deionized water and mobile phase B consisted of 0.1% formic acid in acetonitrile.

Mobile Phase A	Mobile Phase B	Time (min)	Flow Rate (nL/min)
95	5	0	600
10	90	17	600
10	90	29	600
95	5	30	600
95	5	44	600




Figure S1. Proton nuclear magnetic resonance spectrum for TAHI.




Figure S2. Fragmentation spectrum of precursor ion $[M + H]^+$ for TAHI (m/z 427.3) obtained by nano-UPLC-ESI-MS/MS operated in positive ion-mode with nano-electrospray ionization (scan range, m/z 100-500; collision energy, 25 eV).




Figure S3. Fragmentation spectrum of precursor ion $[M + H]^+$ for TAHI (m/z 427.3) obtained by nano-UPLC-ESI-MS/MS operated in positive ion-mode with nano-electrospray ionization (scan range, m/z 100-500; collision energy, 35 eV).




Figure S4. Proton nuclear magnetic resonance spectrum for TAAHI.




Figure S5. Carbon-13 nuclear magnetic resonance spectrum for TAAHI.




Figure S8. Proton nuclear magnetic resonance spectrum for TAHpI.




Figure S9. Fragmentation spectrum of precursor ion $[M + H]^+$ for TAHPi (m/z 469.3) obtained by nano-UPLC-ESI-MS/MS operated in positive ion-mode with nano-electrospray ionization (scan range, m/z 100-500; collision energy, 25 eV).




Figure S10. Fragmentation spectrum of precursor ion $[M + H]^+$ for TAHPi (m/z 469.3) obtained by nano-UPLC-ESI-MS/MS operated in positive ion-mode with nano-electrospray ionization (scan range, m/z 100-500; collision energy, 35 eV).