

## Pfizer-BioNTech COVID-19 Vaccine

**Delivery Checklist** 





COVID-19 vaccine (Pfizer) requires ultra-cold temperatures during shipment and will arrive in a thermal shipping container with dry ice. To unpack the vaccine, staff will need proper personal protective equipment (PPE) and should know how to handle dry ice safely. Training materials from the manufacturer, including videos, are available at www.cvdvaccine.com. Use the following checklist

to ensure the vaccine is received, stored, and handled appropriately.		
1.	Vaccine will be delivered by FedEx or UPS.	8. Wearing PPE, continue to unpack the shipping container following the manufacturer's guidance. Materials
2.	Examine the shipping container for signs of physical damage.	outlining the unpacking process can be found inside the
		container immediately after opening the lid.
	If the thermal shipping container appears	Inspect the tray(s) of vaccine for damage.
	damaged — <b>STOP.</b> Contact Pfizer with questions	Ensure the correct number of trays were delivered.
	or concerns.	Place the tray(s)or box(es) with the vials upright in the
	If the thermal shipping container appears to be	storage unit.
	in good condition, continue to Step 3.	<ul> <li>If storing vaccine at ultra-cold temperatures between -80°C and -60°C (-112°F and -76°F), do NOT open the</li> </ul>
3.	Place the thermal shipping container on the	tray(s) or touch the vials. Return tray to ultra-cold
	floor to unpack (may weigh up to ~80 pounds).	storage within <b>5 minutes</b> .
		<ul> <li>If storing the vaccine at freezer temperatures between</li> </ul>
4.	The thermal shipping container contains dry	-25°C and -15°C (-13°F to 5°F), the vaccine may be
	ice	stored for up to 2 weeks.
	Before opening:	» Vials stored in the freezer may be returned one time to ultra-cold temperature storage (-80°C to -60°C [-112°F
	Collect PPE, including goggles and insulated gloves.	to -76°F]).
	Place the thermal shipping container in an area	Once returned to ultra-cold storage, the 2-week time
	with proper ventilation. A leak in a confined area may create an oxygen-deficient environment and	frame is suspended.
	result in suffocation.	If storing the vaccine at refrigerator temperatures
	result in surfaceation.	between 2°C and 8°C (36°F and 46°F), the vaccine must
<b>5.</b>	Open the thermal shipping container.	be used within 1 month (31 days).
		"> Thawed vaccine cannot be refrozen.
6.	Press and hold the "Stop Shipment" button	Dispose of the dry ice according to the manufacturer's directions.
	on the temperature monitoring device (TMD)	directions.
	for 5 seconds.	9. If using the thermal shipping container to store
	This triggers an e-mail report from the	vaccine,* replenish the container with dry ice pellets
	manufacturer on the temperature status of the	(sized 10 mm to 16 mm) within 24 hours of delivery.
	container during transit. The report will be sent to the provider (facility) e-mail address associated	Unless you have opted out of receiving dry ice when
	with the order.	the order was placed, dry ice will be provided for the
	■ The LED indicator light on the TMD will change	first re-icing.† Close the container using packing tape.
	from blinking to a solid light.	10. Respond to the manufacturer's status report email
	-	regarding ongoing temperature monitoring of the
<b>7.</b>	Proceed based on the color of the	thermal shipping container.
	LED indicator light.	If storing vaccine in an ultra-cold freezer or refrigerator, click
	Green: Vaccine can be unpacked.	on the link to opt out of ongoing temperature monitoring.
	Continue to step 8.	If storing vaccine in the thermal shipping container,* an
	■ <b>Red or no color: </b> STOP! Wait for the status report	additional e-mail will be sent. Add additional contacts to
	on the vaccine.	be notified of the temperature status of the container.
	<ul> <li>Contact the manufacturer immediately if the</li> </ul>	Include after-hours phone numbers.
	status report indicates a temperature excursion	* CDC recommends using the thermal shipping container
	has occurred.	for temporary storage only.

for temporary storage only.

<sup>†</sup> Dry ice will not be provided for Pfizer 450 COVID-19 vaccine orders.

• Continue to step 8 if the vaccine can be used.