

Supplementary Appendix S1—*Leptospira* spp serovars included on MAT panels by diagnostic laboratory.

Diagnostic Laboratory	Autumnalis	Bratislava	Canicola	Grippotyphosa	Hardjo*	Icterohaemorrhagiae	Pomona
CDC Laboratory†	x	x	x	x	x	x	x
Commercial Laboratory A	x	x	x	x		x	x
Commercial Laboratory B	x	x	x	x	x	x	x
University Laboratory A		x	x	x	x	x	x
University Laboratory B			x	x	x	x	x

*Or a serovar in the same serogroup (Sejroe); CDC laboratory uses serovar Wolfi on the MAT panel, which is in the same serogroup as Hardjo.

†Additional serovars tested: *L. interrogans* serovars Australis, Bataviae, Djasiman, Mankarso, Pyrogenes, and Wolfii (Sejroe serogroup, which is the same serogroup as serovar Hardjo), *L. borgpetersenii* serovars Ballum, Javanica and Tarassovi, *L. weilii* Celledoni, *L. kirschneri* Cynopteri, *L. santarosai* serovars Alexi, Borincana and Georgia.