Title

Multi-Walled Carbon Nanotubes Induce Arachidonate 5-Lipoxygenase Expression and Enhance the Polarization and Function of M1 Macrophages - Dataset

Introduction

Fibrogenic carbon nanotubes (CNTs) induce the polarization of M1 and M2 macrophages in mouse lungs. Polarization of the macrophages regulates the production of proinflammatory and pro-resolving lipid mediators (LMs) to mediate acute inflammation and its resolution in a timedependent manner. Here we examined the molecular mechanism by which multi-walled CNTs (MWCNTs) induce M1 polarization in vitro, with a focus on the induction of arachidonate 5lipoxygenase (Alox5), a key enzyme in the biosynthesis of LMs. Treatment of J774A.1 murine macrophages with MWCNTs increased the expression of Alox5 mRNA in a concentration- and time-dependent manner, with the largest induction (5.3-fold over control) occurring at 10 mg/ml and 3 days post-exposure. The Alox5 protein was similarly induced. Alongside Alox5 induction, MWCNTs induced the expression of CD68, a cell surface marker of M1, by approximately 2.0fold and that induction persisted for up to 3 days post-exposure. Both the expression and activity of inducible nitric oxide synthase, an intracellular marker of M1, were increased by MWCNTs up to 8.4-fold. Consistent with M1 polarization, MWCNTs induced the production and secretion of proinflammatory cytokines tumor necrosis factor-α and interleukin-1β, and proinflammatory LMs leukotriene B4 (LTB4) and prostaglandin E2 (PGE2). Moreover, cell-free media from MWCNTpolarized macrophages induced the migration of neutrophils, which was blocked by Acebilustat, a specific leukotriene A4 hydrolase inhibitor, or LY239111, a LTB4 receptor antagonist, but not the cyclooxygenase 2 inhibitor, NS-398, revealing LTB4 as a major mediator of neutrophil chemotaxis from MWCNT-polarized macrophages. Knockdown of Alox5 using specific small hairpin-RNA suppressed MWCNT-induced M1 polarization, LTB4 secretion, and migration of neutrophils, implicating an important role of Alox5 in MWCNT-induced M1 polarization and function. Taken together, these findings demonstrate the polarization of M1 macrophages by MWCNTs in vitro and highlight induction of Alox5 as an important mechanism by which MWCNTs promote proinflammatory responses by boosting M1 polarization and production of pro-inflammatory LMs.

Methods Collection

1. Dispersion of MWCNTs

- MWCNTs (Mitsui-7) were prepared in a control medium (CM) [Dulbecco's modified eagle's medium (DMEM) with 1% fetal bovine serum (FBS)] at a concentration of 2 mg/ml through vortexing and by sonication.
- Stock solutions of MWCNTs were further diluted with the culture media at 2.5 or 10 μg/ml, sonicated immediately before use.

2. Cell culture, polarization, and treatment

- J774A.1 murine monocyte/macrophage cells were grown in DMEM with 10% fetal bovine serum.
- J774A1. Cells (5 x 10⁵ cells/ml) were seeded in DMEM with 3% FBS for 1 day and then treated with other reagents for 1day or 3 days.

- M1 polarization was induced with interferon γ (20 ng/ml) plus lipopolysaccharides (100 ng/ml) and M2 polarization was induced with interleukin 4 (20 ng/ml) for indicated time (typically three days).
- Control media (DMEM media with 1% FBS) were prepared and treated to establish a negative control response.
- After treatment with a single dose of MWCNTs (2.5 or 10 μg/ml), the cells were incubated for 1 day or 3 days in the same media until harvested to examine total RNA and protein expression.
- HL-60, a promyelocytic cell line, was grown in the RPMI1640 medium with 10% FBS.
- Differentiated HL-60 cells (dHL-60) were prepared with use of 2 μM all-trans retinoic acid treatment for 3 days.

3. Quantitative real-time PCR (RT-qPCR)

- Detection and quantification of Alox5, Alox5ap, or β-actin at mRNA level.
- Real-time qPCR was performed, and the relative expression change values were calculated as $2^{-\Delta\Delta Ct}$ and expressed as fold change in comparison with untreated control.

4. Immunoblotting

- Detection and quantification of Alox5, CD68, or β-actin at protein level.
- Images were scanned using HP scanjet and were used to quantify each band with ImageJ software.

5. Detection of nitric oxide synthase 2 (Nos2, mouse iNOS) expression

- Detection and quantification of Nos2 expression using an intracellular Nos2 detection assay kit.
- The fluorescence signal was measured using a fluorescence plate reader and the relative fluorescence unit (RFU) was expressed.

6. Enzyme-linked immunosorbent assay (ELISA)

• Production of proinflammatory cytokines (TNF-α, IL-1β) and lipid mediators (LTB4, PGE2) in cell-free culture supernatant were measured using ELISA kits with a microplate reader equipped with SOFTmax PRO 4.0.

7. Chemotaxis assay

- Transwell cell migration assay was performed using a cell migration assay kit to determine the effect of MWCNTs on neutrophilic cell migration *in vitro*.
- ATRA-induced differentiated HL-60 (dHL-60) cells were collected and suspended in RPMI 1640 medium without FBS (2.5×10⁵ cell/well/100 μl) and plated on each upper chamber.
- The lower chambers were filled with cell-free culture supernatants, indicated amount of MWCNTs in the basal RPMI 1640 medium containing 1% FBS, or RPMI 1640 medium containing 10% FBS.
- For inhibition assay, a cell-free culture media from cells treated with a specific cyclooxygenase 2 inhibitor, NS-398 (at 2 or 10 μM), a leukotriene A4 hydrolase inhibitor, Acebilustat (at 1 or 5 μM), or dimethyl sulfoxide (DMSO) as a vehicle 6 hr prior to MWCNTs or IFN-γ+LPS exposure added at each lower chamber. A specific LTB4 receptor inhibitor LY293111 (5 or 25 nM) or DMSO as a vehicle was added directly into the dHL-60 cell suspension.

• The fluorescence signal was measured using a fluorescence plate reader and the relative fluorescence unit was expressed to show the migrated cells.

8. Alox5 gene silencing

- To knockdown Alox5 gene expression in macrophages, mouse Alox5-specific short hairpin (shRNA) as RNA interference was introduced into macrophages and incubated for 2 days and then cells were used for RT-qPCR and immunoblotting analysis.
- Cells were treated with MWCNTs or M1 inducer and subsequently used for RT-qPCR and immunoblotting analysis, and the cell-free culture media was collected and used for ELISA assay.

Citations: Publications based on the dataset

Lim CS, Veltri B, Kashon M, Ma Q (2022) Multi-Walled Carbon Nanotubes Induce Arachidonate 5-Lipoxygenase Expression and Enhance the Polarization and Function of M1 Macrophages. Nanotoxicology, DOI: 10.1080/17435390.2023.2204161.

Acknowledgements

This work was supported by NIOSH. CAN# 09390BMX.

Authors

Chol Seung Lim¹ (<u>oxr9@cdc.gov</u>)
Brandon Veltri² (<u>otu3@cdc.gov</u>)
Michael Kashon³ (<u>mqk1@cdc.gov</u>)
Qiang Ma^{1,4} (<u>qam1@cdc.gov</u>)

- ¹ Receptor Biology Laboratory, Toxicology and Molecular Biology Branch, HELD, NIOSH, Morgantown, WV 26505, USA
- ² Department of Microbiology, Immunology, and Cell Biology, Robert C. Byrd Health Sciences Center, West Virginia University, Morgantown, WV 26505, USA
- ³ Bioanalytics Branch, HELD, NIOSH, Morgantown, WV 26505, USA
- ⁴ Corresponding author

Contact

For further information, contact: TMBB/HELD/NIOSH 304.285.6241