## Supporting information for

Development of a nanobody-alkaline phosphatase fusion protein and its application in a highly sensitive direct competitive fluorescence enzyme immunoassay for detection of ochratoxin A in cereal

AUTHOR NAMES. Xing  $Liu^{\dagger,\,\,\dagger}$ , Yang  $Xu^{*,\,\,\dagger,\,\,\dagger}$ , De-bin  $Wan^{\S}$ , Yong-hua Xiong $^{\dagger,\,\,\dagger}$ , Zhen-yun  $He^{\dagger,\,\,\dagger}$ , Xian-xian  $Wang^{\dagger,\,\,\dagger}$ , Shirley J.  $Gee^{\S}$ , Dojin  $Ryu^{\sharp}$ , and Bruce D.  $Hammock^{\S}$ 

## **AUTHOR ADDRESS.**

<sup>†</sup>State Key Laboratory of Food Science and Technology, Nanchang University, No. 235 Nanjing East Road, Nanchang 330047, P. R. China

<sup>‡</sup>Sino-Germany Joint Research Institute, Nanchang University, No. 235 Nanjing East Road, Nanchang 330047, P. R. China

§Department of Entomology and UCD Comprehensive Cancer Center, University of California, Davis, California 95616, United States

\*School of Food Science, University of Idaho, 875 Perimeter Drive, MS 2312, Moscow, Idaho 83844-2312, United States




Figure S-1. Colony PCR analysis of recombinant plasmid pecan45-Nb28-AP. Lane M, 100 bp DNA ladder; Lane 1-10, colony PCR products of ten single transformed colonies randomly picked from the Luria-Bertani (LB)-agar plates containing 50 μg/mL carbenicillin. Each PCR reaction (20 μL) consisted of 10 × high fidelity PCR buffer, 10 mM dNTP mixture, 50 mM MgSO<sub>4</sub>, 10 μM of primer AP-F and AP-R, transformed bacteria, platinum Taq high fidelity, and sterilized water. PCR was carried out as follows: 94 °C for 2 min, followed by 30 cycles of 94 °C for 30 s, 55 °C for 30 s, and 68 °C for 1 min. The 433-bp target DNA was confirmed by agarose gel electrophoresis.

Table S-1. Sequence of primer

Name	Sequence (5'-3')
AP-F	CATGCCATGACTGTGGCCCAGCCGGCCCAGKTGCAGCTCG
	TG GAGTCNGGNGG
AP-R	CATGCCATGACTCGCGGCCCCCGAGGCCTCGCCTTGTGGT
	TTTGG TGT CTT GGG

Table S-2. Precursor and corresponding product ions for the MRM detection

Chan Reaction	Dwell (secs)	Cone Volt.	Col. Energy
404.00 > 239.00	0.100	28.0	22.0