



HHS Public Access

Author manuscript

J Biomol Screen. Author manuscript; available in PMC 2017 December 04.

Published in final edited form as:

J Biomol Screen. 2014 February ; 19(2): 325–330. doi:10.1177/1087057113497799.

Profiling Cholinesterase Adduction: A High-Throughput Prioritization Method for Organophosphate Exposure Samples

Melissa D. Carter¹, Brian S. Crow¹, Brooke G. Pantazides², Caroline M. Watson², B. Rey deCastro¹, Jerry D. Thomas¹, Thomas A. Blake¹, and Rudolph C. Johnson^{1,*}

¹Division of Laboratory Sciences, National Center for Environmental Health, Centers for Disease Control and Prevention, Atlanta, GA

²Oak Ridge Institute for Science and Education Fellow at the Centers for Disease Control and Prevention, Atlanta, GA

Abstract

A high-throughput prioritization method was developed for use with a validated confirmatory method detecting organophosphorus nerve agent exposure by immunomagnetic separation-HPLC-MS/MS. A ballistic gradient was incorporated into this analytical method in order to profile unadducted butyrylcholinesterase (BChE) in clinical samples. With Zhang, *et al.* 1999's Z'-factor of 0.88 ± 0.01 (SD) of control analytes and Z-factor of 0.25 ± 0.06 (SD) of serum samples, the assay is rated an "excellent assay" for the synthetic peptide controls used and a "double assay" when used to prioritize clinical samples. Hits, defined as samples containing BChE Ser-198 adducts or no BChE present, were analyzed in a confirmatory method for identification and quantitation of the BChE adduct, if present. The ability to prioritize samples by highest exposure for confirmatory analysis is of particular importance in an exposure to cholinesterase inhibitors such as organophosphorus nerve agents where a large number of clinical samples may be collected. In an initial blind screen, 67 out of 70 samples were accurately identified giving an assay accuracy of 96% and yielded no false negatives. The method is the first to provide a high-throughput prioritization assay for profiling adduction of Ser-198 BChE in clinical samples.

Keywords

butyrylcholinesterase; high-throughput screening; organophosphorus nerve agent; cholinesterase inhibitors; sample prioritization; ballistic gradient

Introduction

Cholinesterases are a family of enzymes comprised of acetylcholinesterase and pseudocholinesterase, also known as butyrylcholinesterase (BChE).^{1,2} These enzymes are

*Correspondence to: R.C. Johnson, Division of Laboratory Sciences, Centers for Disease Control and Prevention, Atlanta, GA 30341, USA. rmj6@cdc.gov.

Disclaimer

The findings and conclusions in this article are those of the authors and do not necessarily represent the views of the Centers for Disease Control and Prevention. Use of trade names is for identification only and does not imply endorsement by the Centers for Disease Control and Prevention, the Public Health Service, or the US Department of Health and Human Services.

responsible for triggering hydrolysis of the neurotransmitter acetylcholine.^{3,4} Hydrolysis of acetylcholine follows nerve stimulation and allows the activated neuron to return to its resting state. Inhibitory compounds like organophosphorous nerve agents can block this process by binding to the cholinesterase serine active site causing a buildup of acetylcholine.^{4,5} This in turn causes the neuron to remain stimulated, sending repeated nerve impulses and leading to constant neuromuscular contraction. In this state, the nervous system is suppressed and patients experience a variety of symptoms related to cholinergic inhibition.⁵

When patients exhibit symptoms of cholinergic inhibition, their cholinesterase activity levels are quickly screened so that treatment, if needed, can be administered immediately. Initial screens of activity are effective for fast diagnoses but lack the specificity needed to identify and confirm the chemical agent of exposure.⁴ Confirmatory methods undergo formal validation and are reviewed annually within a quality assessment system. The methods to detect organophosphate exposure range from GC-MS⁶ and MALDI-TOF-MS⁷ to HPLC-MS/MS.⁷⁻¹¹ For each of these methods, the time-to-first-result can vary due to availability of reference materials, the number of analytes monitored and the complexity of the sample matrix. These limiting factors can be further affected by the number of samples received for analysis.

While confirmatory methods ensure high quality analyses for all samples, little discussion on prioritizing those samples has been available until recently.¹² By profiling a validated method analyte known to be affected by exposure, analysts can organize samples so that those suspected of highest exposure are analyzed first. The following HPLC-MS/MS analysis using a ballistic gradient is designed to complement the validated confirmatory method currently in use. This high-throughput assay is used to profile BChE Ser-198 adduction^{8,10} in serum samples and consequently, provide an effective first tier tool for prioritizing samples for further analysis.

Materials and Methods

Materials

Synthetic unlabeled and stable isotopically-labeled BChE peptides (FGESAGAAS and ¹³C₉-FGESAGAAS, respectively) of 90% purity were obtained from TNO (The Netherlands). BChE monoclonal antibodies from clone 3E8 were commercially available from ThermoFisher Affinity Bioreagents (Rockford, IL, USA). HPLC grade solvents acetonitrile and deionized water were purchased from Tedia (Fairfield, OH, USA). Formic acid (98%), phosphate buffered saline with Tween 20 (PBST) dry powder, dimethyl pimelimidate dihydrochloride (DMP), 0.2 M Tris buffered saline (TBS) 10X concentrate, triethanolamine buffer solution and pepsin from porcine gastric mucosa were all ordered from Sigma Aldrich (St. Louis, MO, USA). Dynabeads protein G were purchased from Life Technologies (Carlsbad, CA, USA). KingFisher 96 Flex microplates (200 µL), KingFisher 96 tip combs for deep well magnets, KingFisher Flex microtiter deepwell 96 plates (v-bottom), protein precipitation plates, 96-well PCR plates (ABgene product), PCR foil and easy pierce 20 µm heat sealing foil were all obtained from Fisher Scientific (Rockford, IL, USA).

Sample preparation

Samples were processed as previously described and recently modified.^{9,10} Briefly, 125 μL of each serum sample was pipetted into a multiscreen HTS HV 0.45 μm opaque filter plate, sealed with adhesive PCR foil, and centrifuged at $3,000 \times g$ for 5 minutes at 20 $^{\circ}\text{C}$ to remove any fibrous tissue. A 75 μL aliquot of the filtered serum was transferred to a 96-well KingFisher deep well plate. A 50 μL aliquot of prepared BChE antibody-coated magnetic beads was placed into all wells of a 96-well KingFisher shallow well plate that corresponded to serum-containing wells. A ThermoScientific KingFisher Flex magnetic particle processor was used to transfer the BChE antibody-coated magnetic beads to the wells containing serum. BChE protein was allowed to bind to the antibody beads by shaking with an Eppendorf MixMate at 1,400 rpm for 2 hours at room temperature. The protein bound beads were transferred using the KingFisher to a plate containing pepsin for protein digestion. The pepsin plate contained 10 μL of 850 ppb internal standard (isotopically labeled synthetic unadducted BChE peptide ($^{13}\text{C}_9$ -FGESAGAAS) and 75 μL 0.6% formic acid. The digestion plate was heated at 37 $^{\circ}\text{C}$ for 30 minutes on an Eppendorf Thermoshaker mixing at 1000 rpm for 10 seconds per minute. The beads were removed from the sample following digestion using the KingFisher Flex. A 285 μL aliquot of acetonitrile was added to a Pierce 2 mL protein precipitation plate, 0.2 μm filter. The samples from the digestion plate were then transferred by manual pipetting. The samples were filtered through the protein precipitation plate using a vacuum manifold, eluting into a 2 mL 96-well Nunc plate. The samples were dried to completeness under nitrogen pressure using a TurboVap at 70 $^{\circ}\text{C}$. The samples were resuspended in 75 μL of 0.6% formic acid and transferred to a 96-well PCR plate. The plate was heat sealed with foil for HPLC-MS/MS analysis.

Preparation of stock solutions

Quality control (QC) pooled human serum was obtained from Battelle Memorial Institute (Columbus, OH, USA) and characterized for use.^{9,10} QCs were used to set the threshold for the prioritization method. The QC materials consisted of unexposed serum and VX exposed serum.

HPLC-MS/MS

BChE levels in human serum were determined on an AB Sciex 6500 triple quadrupole mass spectrometer (AB Sciex, Framingham, MA, USA). The instrument was tuned and calibrated monthly over a mass range of m/z 118–2720 using an Agilent ESI tuning mixture (P/N G2421A) and manufacturer autotune procedure. HPLC elution was performed using an Agilent 1290 Infinity series HPLC system (Agilent, Santa Clara, CA, USA). Samples (2 μL) were injected onto a Waters Acquity UPLC BEH C18 VanGuard PreColumn (2.1 \times 5 mm, 1.7 μm) (Waters, Milford, MA, USA). PreColumn and autosampler temperatures were 25 $^{\circ}\text{C}$ and 10 $^{\circ}\text{C}$, respectively. Mobile phases were 0.1% formic acid in (A) water/acetonitrile (99.9:0.1) and in (B) water/acetonitrile (0.1:99.9). The following ballistic gradient conditions were used at a constant 250 μL flow rate and average back pressure of 143 bar: 0 min 0%B; 0–0.01 min linear gradient 100%B; 0.01–0.35 min 100%B; 0.35–0.36 min linear gradient 0%B; 0.36–0.50 min 0%B. The following optimized instrument parameters were applied for the detection of the unlabeled analyte and the isotopically-labeled internal

standard: collision gas = 12 psig; curtain gas = 10 psig; ion source gas 1 = 40 psig; ion source gas 2 = 60 psig; ion spray voltage = 5500 V; temperature = 500 °C; collision exit potential = 20 V; declustering potential = 136 V; entrance potential = 10 V; dwell time = 35 ms; 5.0 mm vertical probe height and a 'unit' resolution of 0.7 amu at full width half max. Detection of the BChE pepsin digest product peptide FGESAGAAS was determined by selected reaction monitoring (SRM) (BChE detection ion m/z 796.3 \rightarrow 691.3, collision energy = 34 V; confirmation ion m/z 796.3 \rightarrow 620.3, collision energy = 36 V and $^{13}\text{C}_9$ -BChE internal standard m/z 805.4 \rightarrow 700.3, collision energy = 34 V) in ESI positive ion mode.

Data acquisition and processing

Crystal structures were produced using PyMOL v. 1.5 from Schrödinger LLC (Portland, OR, USA). Data acquisition and spectral analysis were carried out utilizing AB Sciex Analyst v. 1.6 build 6536. The capacity factor for the separation was calculated by $k' = (t_r - t_m) / t_m$, where t_r is the analyte retention time and t_m is the system void volume. Zhang, *et al.* 1999's Z -factor of the 96-well plate assay was calculated as $Z = 1 - ((3\sigma_s + 3\sigma_c) / |\mu_s - \mu_c|)$, where σ is the standard deviation, μ is the mean, c is the control, s is the sample and abs is the absolute value.¹³ Similarly, the Z' -factor was calculated.¹³ Percent relative standard deviation $\%RSD = (SD / C_{avg}) \times 100$ was calculated as a measure of assay precision, where C_{avg} is the average concentration calculated, and SD is the standard deviation of C_{avg} . Intraplate precision was determined with an analyte $n = 48$ and interplate accuracy was determined with a plate $n = 3$. Assay accuracy was calculated as $\%Accuracy = (n_{corr} / n_{total}) \times 100$, where n_{corr} are the number of samples prioritized correctly and n_{total} are the number of samples total. The assay's confirmation rate was determined as a percentage by $\%CR = (\text{hits confirmed} / \text{hits selected}) \times 100$, where hits are the samples found to have low unadducted BChE levels and thereby high priority for further analysis.

Convenience sample set

A set of 100 individual human sera samples was commercially obtained from Tennessee Blood Services (Memphis, TN, USA).

Results and Discussion

A novel validated method for extracting human BChE from serum or plasma by immunomagnetic separation coupled with isotope-dilution HPLC-MS/MS was developed by Sporty and colleagues in 2010.^{8,10} This combination of analytical techniques enabled analysts to obtain and quantitatively measure a direct protein biomarker of organophosphate exposure from clinical samples. Since BChE is functional in 99.9% of the general population,¹⁴ profiling this biomarker provides a meaningful target for prioritizing samples received during a potential exposure to cholinesterase inhibitors. Detection of organophosphate adducts to the BChE target are projected up to 55 days post-exposure, expanding response time by 95% from the urinary metabolite method's three day response time.^{6,11} Although further automation of the adduct method led to more efficient sample processing times, sample reporting remained dependent on a 10 hr/96-well plate HPLC-MS/MS analysis and the subsequent data review of the multi-analyte system.⁹ If samples

were prioritized before they were analyzed with the more complex confirmatory methods, samples with the highest exposure profiles would have a faster time-to-report. The benefits of this prioritization design ensure that investigators obtain information about a chemical emergency in the order of highest exposure first.

The key to adding a prioritization method before the confirmatory method was to not prolong the original time-to-first-result. The confirmatory method analysis time for initial calibration and QCs was 1.75 hr. A high-throughput prioritization method was designed to profile at least 96 unknown clinical samples in less than 1.75 hr. The desired method was achieved using a ballistic HPLC gradient. First termed in 2000, ballistic gradients are rapid gradients at high flow rates on short and narrow-bore columns.¹⁵ A 0.5 min ballistic gradient HPLC-MS/MS analysis was developed to target unadducted BChE peptide at a fifth of the confirmatory method's sample volume (Figure 1). A BChE peak at 0.24 min was observed with a measured $k' = 3.8$. Since the BChE peptide used was the same as the biomarker peptide extracted for the confirmatory method, the method had the additional benefit of not altering the original sample processing protocol.¹⁰ The method enabled BChE profiles to be collected from a 96-well plate in 1.6 hr and queued into the confirmatory analysis method in order of suspected exposure.

The analysis of the unadducted BChE in each sample included monitoring the ESI-MS/MS target ion transition m/z 796.3 \rightarrow 691.3; confirmatory target ion transition m/z 796.3 \rightarrow 620.3 and the isotopically-labeled internal standard ion transition m/z 805.4 \rightarrow 700.3. The peak area ratio of the target to the internal standard was used to establish assay parameters. Assay quality was determined using Zhang, *et al.* 1999's high-throughput Z-factor parameter.¹³ A checkerboard 96-well plate of alternating unadducted and adducted synthetic BChE peptide standards were used to assess the developed method's ability to discriminate between unadducted BChE target and adducted BChE (Figure 2A). The peptide standards were prepared at 37 ng/mL, equivalent to the QC material concentration measured in the confirmatory method. Only baseline detection of adducted BChE synthetic standards was observed due to the high specificity of MS/MS analysis. The calculated Z' -factor was 0.88 ± 0.01 (SD) and described an "excellent assay" capable of distinguishing the target with confidence.¹³ The high precision of the peptide standard was 4 %RSD intraplate ($n=48$) and 7 %RSD interplate ($n=3$). A naturally lower Z-factor was predicted and observed for matrix serum targets from clinical samples (Figure 2B). The calculated Z-factor for serum QC reference material was found to be 0.25 ± 0.06 (SD) and termed a "double assay".¹³ This classification was as intended since the profiling method was designed as a means of prioritization before the validated confirmatory method. Profiles of unadducted BChE extracted from the unexposed QC reference material were found to be less precise than the synthetic standards but still within a high working range of 13 %RSD intraplate ($n=48$) and 15 %RSD interplate ($n=3$).

A high confidence in the initial assay quality was established, and a commercial convenience set of individual human serum samples were tested. The BChE target peptide was observed in all of the 96 individual samples with an average response not unlike that observed in the serum QC reference materials (Figure 2C). An intraplate precision of 19 %RSD was calculated among the individual samples and was used to set the threshold for the assay. A

threshold of 19 % of the mean was representative of the reported 22 % population variability.⁵ Samples with a response lower than one standard deviation of the mean (equivalent to 19 %RSD) would be classified as “hits” or high priority samples in the confirmatory analysis queue. The authors understand a threshold at one standard deviation of the mean will return an abnormally high number of false positives due to the inherent target BChE variability in the general population, but this low threshold will also reduce, if not eliminate, the number of false negatives observed. The importance of avoiding false negatives in a chemical emergency or exposure sample set is of utmost importance for formal investigation purposes.

The final prioritization method parameters were applied to a blind response sample set. A blinded plate containing 70 potential exposure samples was prepared and profiled for priority of suspected organophosphate exposure. Samples were comprised of an array of QC materials and individual serum samples. Of the 70 samples analyzed, 67 samples were classified correctly giving a high assay accuracy of 96%. The confirmation rate of the identifications was also high at 93% which means that of the 45 hits prioritized, 42 were later confirmed (Figure 3). The three unconfirmed hits were labeled false positives; no false negatives were observed. A tested assay threshold that did not return false negatives is extremely important when coupled to a confirmatory method for organophosphate exposure.

The described high-throughput prioritization method was developed for use with a validated confirmatory method detecting organophosphate exposure by immunomagnetic separation-HPLC-MS/MS. While all potential exposure samples received would be analyzed with equal care using a validated confirmatory method, sample profiles with suspected organophosphate exposure would be queued for analysis first. A prioritization assay accuracy of 96% and confirmation rate of 93% illustrates the confidence with which this approach can be used. Prioritizing potential exposure samples therefore increases the value of the validated method’s time-to-first-result and will be important in organizing sample analyses in the case of a large scale exposure event.

Acknowledgments

This work was supported by Defense Threat Reduction Agency contract *CDC IAA: OT12-010*. The authors would like to thank Ms. Chariety Sapp of the CDC’s Incident Response Laboratory (IRL) for dispensing convenience set serum samples prior to analysis.

References

1. Çoku ra AN. Butyrylcholinesterase: structure and physiological importance. *Turk J Biochem.* 2003; 28:54–61.
2. Kovarik Z, Radi Z, Berman HA, et al. Acetylcholinesterase active centre and gorge conformations analysed by combinatorial mutations and enantiomeric phosphonates. *Biochem J.* 2003; 373:33–40. [PubMed: 12665427]
3. Valle AM, Radic Z, Rana BK, et al. Naturally occurring variations in the human cholinesterase genes: heritability and association with cardiovascular and metabolic traits. *J Pharmacol Exp Ther.* 2011; 338:125–133. [PubMed: 21493754]
4. Holas O, Musilek K, Pohanka M, et al. The progress in the cholinesterase quantification methods. *Expert Opin Drug Discov.* 2012; 7:1207–1223. [PubMed: 23013366]

5. Nigg HN, Knaak JB. Blood cholinesterases as human biomarkers of organophosphorus pesticide exposure. *Rev Environ Contam Toxicol.* 2000; 163:29–111. [PubMed: 10771584]
6. van der Schans MJ, Polhuijs M, van Dijk C, et al. Retrospective detection of exposure to nerve agents: analysis of phosphofluoridates originating from fluoride-induced reactivation of phosphorylated BuChE. *Arch Toxicol.* 2004; 78:508–524. [PubMed: 15170525]
7. Carletti E, Li H, Li B, et al. Aging of cholinesterases phosphorylated by tabun proceeds through O-dealkylation. *J Amer Chem Soc.* 2008; 130:16011–16020. [PubMed: 18975951]
8. Noort D, Fidder A, Van der Schans M, et al. Verification of exposure to organophosphates: generic mass spectrometric method for detection of human butyrylcholinesterase adducts. *Anal Chem.* 2006; 78:6640–6644. [PubMed: 16970345]
9. Knaack JS, Zhou Y, Abney CW, et al. A high-throughput diagnostic method for measuring human exposure to organophosphorus nerve agents. *Anal Chem.* 2012; 84:9470–9477. [PubMed: 23083472]
10. Sporty JLS, Lemire SW, Jakubowski EM, et al. Immunomagnetic separation and quantification of butyrylcholinesterase nerve agent adducts in human serum. *Anal Chem.* 2010; 82:6593–6600. [PubMed: 20617824]
11. Swaim LL, Johnson RC, Zhou Y, et al. Quantification of organophosphorus nerve agent metabolites using a reduced-volume, high-throughput sample processing format and liquid chromatography-tandem mass spectrometry. *J Anal Toxicol.* 2008; 32:774–777. [PubMed: 19021934]
12. Hamers T, Legler J, Blaha L, et al. Expert Opinion on Toxicity Profiling-Report from a NORMAN Expert Group Meeting. *Integr Environ Assess Manag.* 2013; 9:185–191. [PubMed: 23307398]
13. Zhang JH, Chung TDY, Oldenburg KR. A simple statistical parameter for use in evaluation and validation of high-throughput screening assays. *J Biomolec Screen.* 1999; 4:67–73.
14. Li B, Duysen EG, Saunders TL, et al. Production of the butyrylcholinesterase knockout mouse. *J Molec Neurosci.* 2006; 30:193–195. [PubMed: 17192674]
15. De Nardi C, Bonelli F. Moving from fast to ballistic gradient in liquid chromatography/tandem mass spectrometry pharmaceutical bioanalysis: matrix effect and chromatographic evaluations. *Rapid Commun Mass Spectrom.* 2006; 20:2709–2716. [PubMed: 16921566]




Figure 1. Human BChE was extracted by immunomagnetic separation, digested by pepsin and analyzed by HPLC-MS/MS analysis

The X-ray crystal structure of human BChE with waters removed for clarity was created using PyMOL (PDB identifier 1P0M). The Ser-198 containing BChE peptide (FGESAGAAS) is highlighted in black. A representative extracted ion chromatogram of the peptide transition m/z 796.3 \rightarrow 691.3 is shown at retention time 0.24 minutes ($k' = 3.8$).




Figure 2. Zhang, *et al.* 1999's Z-factor was used to evaluate assay quality across each 96-well plate

The Z-factor parameter was calculated using a checkerboard 96-well plate analysis for the system's **A**) synthetic positive and negative controls (◆-unadducted FGESAGAAS peptide and ▲-VX-adducted FGESAGAAS peptide) yielding an excellent assay Z'-factor of 0.88 ± 0.01 (SD) and **B**) quality control reference materials (◆-human serum and ▲-fully VX-adducted human serum) returning a double assay Z-factor of 0.25 ± 0.06 (SD). **C**) A convenience set of 96 individual serum samples (◆) were analyzed for comparison to the parameters of the quality control materials. The mean (*solid line*) and three standard deviations of the mean (*dashed line*) were calculated for all samples. One standard deviation of the mean (*dotted line*) was set as the assay threshold. Below the threshold, samples were considered high priority or "hits".




Figure 3. The assay was used to prioritize samples from a blind response exercise

The assay was shown to have an accuracy of 96%, evaluating 67 of the 70 samples as expected. There were 42 known hits of high priority (▲) in the sample set; the 42 hits were confirmed plus three false positives (x) and no false negatives, giving a high confirmation rate of 93%. All other samples (◆) are distributed around the mean. The mean (*solid line*), three standard deviations of the mean (*dashed line*) and one standard deviation of the mean (*dotted line*) were the parameters established from the individual serum convenience set.